

SUTURE TRAINING PADS

Mult-Layer Suture Training Pads

The Mult-Layer Suture Training Pads consist of four (4) layers, intended to represent the dermis, subcutaneous layer, linea, and transverse fascia.

Inconsistencies in thickness are derived from the hand crafting process and provide a variation in the training / suturing surface.

Keep pads in bag when not in use.

The pads contain silicone oils which may leach out over time and can stain porous surfaces. Once removed from the bag, the pads should be washed with soap and water to remove the oily residue, if needed. Silicone oil residue can be cleaned off surfaces with a phosphate based detergent and water mixture.

Avoid contact with skin, eyes and porous surfaces.

Hollow Organ Suture Pads

These pads consist of one (1) layer, intended to represent a hollow organ. Inconsistencies in thickness are derived from the hand crafting process and provide a variation in

Canine Suture Pads

These suture pads are the same as our regular multi-layer suture pads, but also feature nipples and umbilicus, which are used as landmarks on the VSI Canine Spay Simulator.

Suture Pad Training Base

The base is supplied with two foam inserts. The solid foam insert is for use with the multilayer suture pad. The perimeter foam insert is for use with the hollow organ suture pad, and can be stored inside the base when not in use.

The suction pads on the bottom of the base will adhere most effectively to a smooth work surface.

The base unit and foam inserts should be cleaned with soap and water. Harsh chemicals or solvents will destroy the plastic and/or foam.

Install the suture training pads by pulling the holes in the pads over the pins in the base.

Please contact Veterinary Simulator Industries Ltd. with specific concerns or inquiries.

*Veterinary Simulator Industries
Calgary, Canada, (403)-262-9393,
www.vetsimulators.com.
consult@vetsimulators.com*

Stay up to date with new products and features on Facebook, YouTube and Instagram [@vetsimulators](https://www.instagram.com/vetsimulators)